

SLIMLINE DESIGN,
LIGHT IN WEIGHT

320
compact



Incorporating the latest technology from Electric Mobility, the Rascal 320 Compact represents a total innovation in lightweight powerchair design. The Rascal 320 Compact is half the size of most powerchairs with the heaviest part weighing only 58lb, making it so easy to take wherever you are going. Quick to dismantle and with rear wheel drive for superb traction, the slimline Rascal 320 Compact will help you live life to the full.

The Rascal 320 Compact gives you –

- convenient offboard charging with a quick release battery pack
- increased independence with a design that is simple to dismantle for easy transportation
- outstanding manoeuvrability indoors and out with a turning circle of only 66cm (26")

- optional extras
lightweight padded seat, elevating legrests, swing-away legrests or dual/attendant control

Combining the benefits of a compact, lightweight design with all the quality and comfort that the Rascal range is renowned for, the Rascal 320 Compact from Electric Mobility is causing a mini powerchair revolution.



POWER WITH PERFORMANCE

Technical Specifications	320 Compact rear wheel drive
Overall length*	67.3cm (26.5")
Overall width	56cm (22")
Maximum rider weight**	136kg (300lb)
Batteries	2x22ah
Battery pack weight	15.2kg (33.5lb)
Total weight	53.3kg (117.5lb)
Heaviest component	26.3kg (58lb)
Maximum range†	20.1km (12.5 miles)
Maximum speed	6.4kph (4mph)
Castor wheel size	17.8cm (7")
Drive wheel size	22.9cm (9")
Turning radius	66cm (26")
Safe climbing angle	8°
Colour	blue
Optional extras	n/a
	Elevating legrests
	Swing-away legrests
	Lightweight padded seat
	Dual/attendant control

Kalabar Mobility (for your mobility needs)
showroom: 37 High Street, Romford, Essex, RM1 1JL
Tel: 01708 755 700

