

Escape PWR SE Powerchair



The Days Healthcare Escape PWR SE Powerchair is a lightweight, 4mph model that is both stylish and rugged. It features compact styling and user-friendly controls making it suitable and effortless to use both indoors and out.

The Days Healthcare Escape PWR SE electric wheelchair is ideal for users looking for independence and portability. You can go shopping, visit friends, or simply go about your daily routine with ease.

The Days Healthcare Escape PWR SE electric wheelchair can be disassembled in seconds so it can be loaded into the boot of your car with ease. You can take it with you wherever you go.

There is a Penny and Giles controller fitted to the Days Healthcare Escape PWR SE powerchair adding to the reliability and making it easy to use.

The Days Healthcare Escape PWR SE electric wheelchair has puncture-proof tyres, alloy rear wheels, and height-adjustable anti-tip wheels, all fitted as standard.

The Days Healthcare Escape PWR SE powerchair has detachable armrests to make side transfers easier. The detachable swing in/out leg rests allow the user to get closer to objects, and to help with front transfers.

The Days Healthcare Escape PWR SE powerchair has a maximum user weight of 18 stone /115Kg and a maximum range of 13 miles.

- Penny and Giles VSI controller
- Puncture-proof tyres
- Alloy rear wheel
- Height adjustable anti-tippers provide safety
- Step-less height adjustable footplate as standard
- Cushioned upholstery

- Premium Quality materials
- Easily folded for transport and storage
- Detachable / Swing-back armrests for ease of transfer
- Detachable Swing-in / Swing-out leg rests allow the user to get closer to objects and aids front transfers
- Heel strap as standard for use comfort and safety
- Manual freewheel function
- Dual and attendant control options
- Elevating legrests (option)

Specifications

Battery charger size 240V - 24V 4Amp amps

Battery size 30A/h Deep-Cycle 12V x 2

Ground clearance 10 cm

Maximum range 13 miles

Maximum speed 4 mph

Safe climbing angle 12 °

Seat depth 41 cm

Seat width 41 cm

Turning radius 95 cm

Weight with batteries 65 kg