



A move in the right direction!

Technology and performance you can count on



Rascal 330 Compass



Technical Specifications

Overall length (inc standard footboard)	97cm (38")
Overall width	62cm (24.5")
Maximum user weight	136kg (300lbs)
Batteries	34 amp
Total weight	83kg (184lbs)
Maximum range	19km (12 miles)
Maximum speed	6.4kph (4mph)
Castor wheel size	15x5cm (6x2")
Drive wheel size	25x9cm (10x3.5")
Turning radius	48cm (19")
Safe climbing angle	6°
Options	50 amp batteries Footboard extension Oxygen bottle holder Swing away legrests Elevating legrests

* Range is calculated under test conditions. Weight of user, terrain and weather conditions can affect range of travel
 ** 50 amp batteries require additional height battery box cover

For further details on the Rascal 330 Compass or to discuss other models in the extensive Electric Mobility range call:

Freephone 0800 252614

Your Electric Mobility dealer:

Kalabar Mobility (for your mobility needs)
 showroom: 37 High Street, Romford, Essex, RM1 1JL
 Tel: 01708 755 700

Due to a policy of continual improvement, Electric Mobility Euro Limited reserves the right to change specifications without prior notice.



Electric Mobility Euro Limited Canal Way Ilminster Somerset TA19 9DL
 T: 01460 258100 F: 01460 258125 E: sales@electricmobility.co.uk www.electricmobility.co.uk

Location photography by kind permission of Somerset College of Arts and Technology



10/12/2006