

## ultralite 760/765



The UltraLite 760 and 765 from Electric Mobility give you the lifestyle you desire. Whichever you choose, you will benefit from a compact 96cm (38") length and a lightweight construction, allowing you to take your powerchair on public transport or in the car without folding the rear seats. Choose between the 760 with its quality padded seat or the 765 which boasts a luxurious, reclining seat for ultimate comfort.

When you want total freedom – choose the UltraLite Powerchair.



The 760 has a padded lightweight seat



The 765 has a luxurious padded seat



The 765's luxurious seat can be reclined and has an adjustable head rest



Handy storage pouch



Folding arm rests aid access to chair



Folds easily for storage



Fits into the boot of most cars



Easy to clean foot plate with carrying handle



Detachable battery pack for easy charging



One handed operation



8" solid tyres with rear wheel drive



Minimal strength required to release the battery pack and seat



Minimal strength required to release the battery pack and seat



8" solid tyres with rear wheel drive



Anti-tip castors for added safety

### Technical Specifications – 760/765

Overall length (folded)	66cm (26")
Overall width	57cm (22")
Maximum rider weight	91kg (200lb)
Batteries	2x12amp
Battery pack weight	10kg (22lb)
Heaviest component weight	19.5kg (43lb)
Seat weight with arms	760 – 10kg (23lb) 765 – 15kg (34lb)
Total weight	760 – 40kg (88lb) 765 – 45kg (99lb)
Maximum range*	14km (9mile)
Maximum speed	6kph (3.8mph)
Front wheel size	13cm (5")
Rear wheel size	20cm (8")
Turning radius	77cm (30")
Safe climbing angle	6°
Colour	Ruby red

Your local Electric Mobility dealer  
**Kalabar Mobility (for your mobility needs)**  
 showroom: 37 High Street, Romford, RM1 1JL  
 Tel: 01708 755 700

For further details on the UltraLite 760 and 765 Powerchairs, or to discuss other models in the extensive Electric Mobility range –

freephone 0800 252614



EASY TO DRIVE, SIMPLE TO FOLD



Electric Mobility Euro Limited, Canal Way, Ilminster, Somerset TA19 9DL  
 Tel: 01460 258100 Fax: 01460 258125 email: sales@electricmobility.co.uk www.electricmobility.co.uk



Due to a policy of continual improvement, Electric Mobility Euro Limited reserve the right to change specifications without prior notice.

07/04