



The Elite Range

For a better quality of living



Features

The level of comfort offered by this Elite riser-recliner becomes even more important for those with a bariatric condition. There's an Elite model to suit every shape and size – the superb 35 stone capacity provides those all essential levels of support and comfort. The three-way action enables you to achieve a truly comfortable position for watching television before reclining the back. And when it comes to standing, the Elite provides you with all the support you need to rise from the chair.



Leg rest extension (optional)



Arm rest covers (standard)



Luxurious head rest pillow (optional)



Diamond Mink



Merlot water-resistant

All Elite models feature:

- Smart Tek
- Arm rest covers
- Removable back rest
- Quick-release hand control
- Magazine pocket
- Battery back-up (power failure)

Also available with:

- Heat and massage
- Leg rest extension
- Head rest pillow
- Castors (lockable)



Elite riser-recliner shown here in Merlot water-resistant fabric

NB Leg-lift capacity 68kg/150lbs

Specifications

	Overall width	Overall height	Seat height	Back rest height	Seat width	Seat depth	Weight limit
Medium	95 cm (37.5")	105.5 cm (41.5")	48 cm (19")	68.5 cm (27")	66 cm (26")	51 cm (20")	227 kg (500 lb)
Large	95 cm (37.5")	109 cm (43")	53.5 cm (21")	71 cm (28")	66 cm (26")	53.5 cm (21")	227 kg (500 lb)
Tall	100 cm (39.5")	117 cm (46")	53.5 cm (21")	79 cm (31")	71 cm (28")	61 cm (24")	227 kg (500 lb)

For further details on the Elite Range or to discuss other models in the extensive Cosi Chair range call:

Freephone 0800 252614

or visit our website:

www.electricmobility.co.uk



Electric Mobility Euro Limited Canal Way Ilminster Somerset TA19 9DL
T: 01460 258100 F: 01460 258125 E: sales@electricmobility.co.uk

Your Cosi Chair dealer:

Due to a policy of continual improvement, Electric Mobility Euro Limited reserves the right to change specifications without prior notice

