

Features

The Ambassador Range is available in three sizes – small, medium and large. Designed and built by specialists in the manufacture of riser-recliners. This outstanding riser-recliner is available in a choice of two fabrics – Diamond (five colours) and Brisa (three colours).



Arm rest covers (standard)



Luxurious head rest pillow (optional)



Leg rest extension (optional)

All Ambassador models feature:

- MaxiComfort™ technology
- Smart Tek
- Arm rest covers
- Removable back rest
- Quick-release hand control
- Magazine pocket
- Battery back-up (power failure)

Also available with:


- Heat and massage
- Leg rest extension
- Head rest pillow
- Castors (lockable)



Brisa Buckskin




Brisa Black Onyx



Brisa Hollyhock




Diamond Sapphire Blue



Diamond Cherry



Diamond Mink



Diamond Paprika



Diamond Jade

Specifications

	Overall width	Overall height	Seat height	Back rest height	Seat width	Seat depth	Weight limit
Small	82.5 cm (32.5")	104 cm (41")	45 cm (17.75")	61 cm (24")	49 cm (19.25")	45 cm (17.5")	136 kg (300 lb)
Medium	82.5 cm (32.5")	109 cm (43")	47 cm (18.5")	65 cm (25.5")	51.5 cm (20.25")	51 cm (20")	170 kg (375 lb)
Large	90 cm (35.5")	114 cm (45")	49.5 cm (19.5")	71 cm (28")	56 cm (22")	54 cm (21.25")	170 kg (375 lb)

For further details on the Ambassador Range or to discuss other models in the extensive Cosi Chair range call:

Freephone 0800 252614

or visit our website:

www.electrictomobility.co.uk



Electric Mobility Euro Limited Canal Way Ilminster Somerset TA19 9DL
T: 01460 258100 F: 01460 258125 E: sales@electrictomobility.co.uk

Your Cosi Chair dealer:

Kalabar Mobility (for your mobility needs)
showroom: 37 High Street, Romford, Essex, RM1 1JL
Tel: 01708 755 700

Due to a policy of continual improvement, Electric Mobility Euro Limited reserves the right to change specifications without prior notice

