

Medicare Enigma XS Standard Aluminium Wheelchair

A high specification wheelchair at an outstanding price, the Enigma XS standard aluminium wheelchair offers an array of added features with a value-for-money price tag which will be hard to beat.

- 18" seat width.
- Lightweight aluminium frame in metallic blue.
- 24" self-propel wheels with quick release mechanism.
- Deep comfort tension adjustable back AND seat canvasses.
- Comes with deluxe seat and back cushions.
- Quick release swing in / swing out leg rest with full height and footplate angle adjustment.
- Flip back and removable armrests.
- Half folding back mechanism.
- Height adjustable anti-tip wheels.
- Adjustable rear wheel position.
- Weighs just 8.3kg / 18lb in carry format, and 16.8kg / 37lb fully assembled.
- 18 stone / 115kg weight limit.

Specification

- Dimension Assembled: (L x W x H) 40" x 25" x 37" / 103cm x 63cm x 93cm.
- Dimension Travel Mode*: (L x W x H) 33" x 10" x 25" / 83cm x 26cm x 63cm.
- Seat (W x D x H): 18" x 17" x 20.5" / 45cm x 42cm x 52cm.
- Weight Assembled: 16.8kg / 37lb.
- Weight Travel Mode*: 8.3kg / 18lb.
- Max. User Weight: 115kg / 18st.

* Refers to the wheelchair folded and detachable parts removed

Wheelchair Accessories



Wheelchair Rucksack
(H13" x W12" x D6")



Wheelchair Storage Bag

As a part of an ongoing product improvement initiative, the manufacturer reserves the right to change specifications and design without notice.