

The Excel Xeryus Position



The Excel Xeryus Position is a unique foldable comfort wheelchair fitted with pressure relieving positioning system, aluminium lightweight wheelchair with foldable frame and combined semi tilt-in-space and recline seat system.

Multi Adjustable wheelchair that takes the worry out off getting the right fit, the optional Position lever (normally fitted under the right armrest) can be used by the patient to adjust his sitting position when a care is not available.

Optional Extras

Headrest: flexible, height adjustable, depth and angle adjustable, foldable

Extra 'Position' adjusting lever: to be operated by user while sitting in the chair

Stump board: detachable swing away, angle adjustable (please state seat width for correct fit)

Features

- Seat height adjustment front, standard 50 cm (3 positions possible)
- Wheelbase adjustment through rear axle plate (4 positions)
- Fitted with pressure relieve positioning system
- Self propel with design rear wheels 24" x 1", PU tyres (Air Tyres optional, no charge)
- Front castors 5" x 1 ¼"(polyurethane) with suspension on front forks, PU Flat free wheels
- Quick-release axles on rear wheels
- Seat fitted with straps, for adjusting seat width and seat shape
- Backrest equipped with 8 tension straps
- Seat and backrest upholstery with pressure relieve material inside
- Bearing house angle adjustable
- Special Order Only

- Removable long armrests, height adjustable, also multiple depth adjustable
- Removable long armrests, height adjustable, also multiple depth adjustable
- Nylon heel loops on foot plates
- Adjustable angle footplate from 90-105 degrees
- Footplates, depth adjustable
- Tension bar with fold down mechanism
- Anti-tip wheels, left/right, height adjustable, removable
- Footrests, swing away, detachable
- Height and angle adjustable push handles
- Cross bar width adjustable, from 40 cm to 50 cm in 2,5 cm steps
- Height adjustable backrest, 2 heights available

Seat width	40, 42½, 45, 47½ and 50 cm
Seat depth	49 cm (can be reduced to 43 cm)
Seat height	48½, 50 and 51½ cm
Backrest height	44 and 46½ cm (adjustable)
Total length	105 cm
Total height	95 cm
Minimal width	33 cm
Armrest height	19-35 cm, height adjustable
Weight (total)	18 kg
Weight largest part	10,2 kg
Delivery Time	2 - 3 weeks
Maximum user weight	130 kg