

DAYS ESCAPE LITE TRANSIT WHEELCHAIR

The Days Medical Escape lite transit wheelchair has been designed to be compact, lightweight, and practical making it ideal for both indoor and outdoor use. Its compact dimensions allow good manoeuvrability indoors, especially when used in confined spaces. With its ultra lightweight construction, the Escape Lite is easy to fold for storage or transportation.

Features

- **Ultra Lightweight**
- **Folding Frame & Backrest**
- **Full Length Armrest Attendant Cable Operated Brakes**
- **Tool Kit/Seat belt Included**
- **Seat Width Option Available in 16 or 18 inch**
- **Swing Away and Detachable Footplates**
- **Available in silver frame only**

Specifications

- **Seat Depth: 41cm**
- **Seat Height: 50cm**
- **Backrest Height: 41cm**
- **Backrest Height Folded: 26cm**
- **Footrest Length: 37cm - 53cm**
- **Armrest Height: 21cm**
- **Overall Length: 90cm**
- **Overall Height: 92cm**
- **Max User Weight: 100kg**
- **Overall Weight: 10kgs**