

Features

- Aluminium lightweight wheelchair with foldable frame
- Under seat pocket (for wallet not all models)
- Crutch Holder (optional no up charge)
- Self propel wit large rear wheels 24 x 1 3/8"
- Flip-backward long armrests (standard)
- Rear wheels, polyurethane (PU), flat free
- Nylon heel straps on foot plates
- Aluminium lightweight wheels with aluminium rear rims
- Adjustable angle footplate from 90-105 degrees
- Front castors 7 x 1 1/4"(polyurethane)
- Tension bar with fold down mechanism
- Quick-release axles on rear wheels
- Height and angle adjustable push handles
- Foldable back and seat upholstery (the tension bar is available on fixed backrests only)
- Tension adjustable seat upholstery
- Backrest equipped with 8 tension straps
- Universal seat cushion rear wheels, 12 1/2" pneumatic and plastic rim
- Seat height adjustment front, standard 50 cm(5 positions possible,of which 3 adjustable and 2 optional at up-charge)
- One seat height front possible, 50 cm
- Wheelbase adjustment through rear axle plate(3 positions, of which 1 amputee)
- Ballhead angle adjustment, also horizontal adjustable
- Depth and angle adjustable aluminium brakes
- Product conforms to standards EN12183
- Product conforms to crash test ISO 7176-19:2001

Weight largest part (folded)	8.4 kg
Maximum user weight	130 kg
Seat width	35, 37 1/2 , 40, 42 1/2, 45 and 50 cm
Seat depth	43 and 46 cm
Seat height	45, 47, 50, 52 and 54 cm
Backrest height	43 cm
Total length	103 cm
Total height	95 cm
Rear wheels	24" x 1 3/8" Rear wheels PU flat free, spoked
Minimal width	31 cm
Models	transit or self propelled
Detachable swing away armrest	yes
Armrest height	26 or (optional) 19-35 cm, height adjustable
Weight (total)	13.9 kg, (transit version: 13.5 kg)